



CREATE/SELECT CONSULTATION




SELECT IMAGES




INPUT RING SIZE & STRUT LENGTH




DETECT FIDUCIALS




CONFIRM THE FRAME



MAP BONE DEFORMITY



CONFIGURE WAYPOINTS




REVIEW THE PRESCRIPTION

PRIMARY IMAGE



SECONDARY IMAGE



FRAME DATA


CONFIGURE DEVICE	
✓ SUPERIOR RING	10160 160mm Full Ring
✓ STRUT 1, mm	251
✓ STRUT 2, mm	275
✓ STRUT 3, mm	224
✓ STRUT 4, mm	230
✓ STRUT 5, mm	201
✓ STRUT 6, mm	254
✓ INFERIOR RING	10160 160mm Full Ring

2D IMAGE DATA



Manual Measurements Unnecessary

Input Known Variables

3D VIRTUAL WORLD




PERSPECTIVE VIEWS



Generates Patient Prescription

Simple • Fast • Easy



AMDT Holdings, Inc.
328 Poplar View Lane East
Suite 2
Collierville, TN 38017

info@amdtholdings.com
(901) 853-4366